
DP-8330

Examination Table







Backrest adjustment (angle: 0-50°)



Depression foot pedal, whole bed can be moved.



Lift up and down (57-86cm)



Adjust the foot pad to keep the bed steady.

PATENT: ZL201520111798.7

MAX LOAD 200KG

strong and stable


Metal frame with laser cutting, robot welding, parts integrating and the overall solid.



This table is simple, good flexibility, safe and stable operation.



Japan imported, bearing 200kg, good sealing performance, no leakage, smooth and good use



Has weak bearing capacity, unstable operation, easy to fail, easy to leak oil and no safety



