MPS/P1

Basic paraffin embedding center

The space saving MPS/P1 is a complete paraffin dispensing unit including a 3.8 liters paraffin container, a heated working area and a cooling spot. The ergonomically elevated heated working area is illuminated by an adjustable LED light and allows comfortable processing of specimen. All functions of the unit are controlled by an intuitive control panel including programmable working hours. Two heated trays provide a pre-warming capacity of 80 cassettes and 250 embedding molds.

Highlights

■ Manual lever and foot switch for paraffin dispensing

- Spacious cooling spot
- Adjustable magnifying glass
- Armrest for optimized work ergonomics
- Spacious and heated paraffin collection tray
- 3.8 liters paraffin container
- Programmable operating times
- Easily adjustable forceps holders
- Adjustable LED illumination
- Intuitive control panel
- Integrated heated wax trimmer:
 - $\quad \hbox{Efficiently removes excess paraffin from any embedding tissue cassette} \\$
 - Paraffin is easily melted away and drains into the removable tray below
 - Physical damages (e.g. due to conventional trimming tools such as knife or razor) are consequently avoided
 - Large scale design of the trimmer area enables simultaneous paraffin trimming of multiple blocks
- Two pre-warmed areas for 80 cassettes and 250 embedding molds.
- Huge variety of electrically heated forceps (optional)

Flexible

Forceps can be freely placed in the holders with a capacity of three forceps each. Electrical forceps can be connected optionally.



Everything within reach

Two heated trays, for either 80 standard cassettes and 250 embedding molds are conveniently located in the back of the working area.



Standard Components

Pre-warming and dispensing module

Foot switch for paraffin dispenser

Magnifying glass, adjustable

Paraffin collection tray

Forceps holder

2 x forceps (1 x blunt, 1 x pointed)

Cleaning spatula black

Base molds 7 x 7 x 7mm, 12 pcs., stainless steel

Base molds 24 x 24 x 7mm, 12 pcs., stainless steel

Mains cable

Operation manual

Specifications

Up to 80°C programmable Temperature range **Cooling spot** 50 x 50mm Paraffin capacity 3.81 Up to 80 cassettes Heated cassette area Heated mold area Up to 250 molds (width x depth x height) 125 x 205 x 40mm each (Tray inserts are removable) Dimensions 450 x 600 x 295mm (width x depth x height) Working area (width x depth) 450 x 240mm

MPS/P1 220–230V, 50/60Hz MPS/P1 115V, 60Hz

Weight without accessories

Item No. 10190000 Item No. 10192000



22ka